

## Making proteins accessible™

Is obtaining soluble, stable and active protein a bottleneck in your research and your next breakthrough?


### Explore eProtein Discovery™

Nuclera empowers scientists to make progress on protein projects through a rapid protein prototyping system to automate construct screening, expression and purification characterization of proteins.

- > **Rapid protein prototyping** enables progress by allowing scientists to gain awareness quickly about which proteins – and which variations of a protein – will work
- > **Simultaneously screen multiple constructs** and protein synthesis reagents for soluble expression, and then scale up to micrograms of recombinant protein off-cartridge to test in your applications
- > Explore multiple DNA constructs, including **solubility tags, truncations, polymorphisms and isoforms** on the same smart cartridge to expand your range of accessible proteins



### eProtein Discovery™ Workflow



## Robust screening data: Soluble expression and purification

Robust solubility screening and purifiable yield assessment provided, allowing for the selection of the best construct and Cell-Free Blend to obtain desired protein.



Nuclera's technology represents fresh approaches, which will improve cost and quality significantly.




### Prof. George Church

Professor of Genetics at Harvard Medical School and Professor of Health Sciences and Technology at Harvard and the Massachusetts Institute of Technology (MIT)


**Figure 1. Expression and purification characterization.** The instrument reports on expression and purified yield (mg/mL or µM) to inform on the most favorable construct and cell-free blend.

## Purified active protein



Lane	Description
1	Molecular weight ladder
2	Cell-free Blend-3 negative control
3	Total protein after cell-free protein synthesis (CFPS)
4	CPFS (soluble) obtained following centrifugation and supernatant retained
5	Unbound; flow through sample not bound to beads
6	Purified (elute); sample eluted off of the beads

**Figure 2. Scale-up expression and purification.** SDS-PAGE showing expressed and purified proteins from a scaled-up reaction.



### Figure 3. Activity testing of VEGF protein produced with eProtein Discovery.

The cell-based PathHunter® dimerization assay was used to observe a change in substrate presence, which was reported in relative light units (RLU) indicating the presence of active protein. Two biological replicates, shown here as a mean, were carried out on two separate occasions. The VEGF protein displayed an EC<sub>50</sub> of 12.49 ng/mL. PathHunter® is a registered trademark of Eurofins DiscoverX as used in US and/or in other countries.

**Key:** ● SUMO\_VEGF\_STREP\_DET ● Negative Control

## Which proteins have been produced so far?



**Figure 4. Proteins produced.** Chaperones, Hydrolases, Ligase, Oxidoreductase, Signaling protein, Structural protein and Transferases with the molecular weight range: Min: 18 kDa to Max: 300kDa (Avg: 46kDa).

nuclera

Cambridge, UK  
Boston, MA USA  
Email: [hello@nuclera.com](mailto:hello@nuclera.com)



[www.nuclera.com](http://www.nuclera.com)